

Research paper

Reliability assessment in advanced nanocomposite materials for orthopaedic applications

Jérôme Chevalier^a, Paola Taddei^b, Laurent Gremillard^a, Sylvain Deville^a, Gilbert Fantozzi^a, J.F. Bartolomé^c, C. Pecharroman^c, J.S. Moya^c, L.A. Diaz^d, R. Torrecillas^d, Saverio Affatato^{e,*}

ABSTRACT

^a Université de Lyon, INSA Lyon, MATEIS - UMR CNRS 5510, 69621 Villeurbanne Cedex, France

^b Centro di Studio sulla Spettroscopia Raman, Dipartimento di Biochimica "G. Moruzzi", Sezione di Chimica e Propedeutica Biochimica,

nanocomposites in the orthopedic field.

Via Belmeloro 8/2, Università di Bologna, 40126 Bologna, Italy

^c ICMM, Spanish Research Council (CSIC), Campus Canto Blanco, Madrid, Spain

^d INCAR, Spanish Research Council (CSIC), La Corredoria s/n Ap. 73, 33080 Oviedo, Spain

^e Laboratorio di Tecnologia Medica, Istituto Ortopedico Rizzoli, Bologna, Italy

ARTICLE INFO

Article history: Received 29 June 2010 Received in revised form 22 October 2010 Accepted 26 October 2010 Published online 2 November 2010

Keywords: Ceramic nano-composites Crack resistance Aging resistance SEM AFM Photoluminescence

1. Introduction

Nowadays, the main issue for THR is the generation of wear debris produced mainly by the acetabular component (Campbell et al., 2004). Standard artificial hip joints consist of a polymer cup made of Ultrahigh Molecular Weight Polyethylene (UHMWPE), placed in the acetabulum via a metal-back component, and a metal (stainless steel or cobalt chromium alloy) ball fixed to a metal stem introduced in the femur. Any use of the joint results in cyclic stress of the polymer cup against the metal ball. During the reciprocating motion of normal joint use, UHMWPE fibrils are then

© 2010 Elsevier Ltd. All rights reserved.

Alumina-zirconia nano-composites were recently developed as alternative bearing ma-

terials for orthopedics. Previous, preliminary reports show that such alumina-zirconia

nanocomposites exhibit high crack resistance and low wear rate. In this paper, additional

information is given in terms of wear, crack resistance and ageing behaviour: femoral heads

are inspected after 7 million cycles of wear testing on a hip simulator, crack resistance is measured and compared to other ceramics used today in orthopedics, slow crack growth

is reported under static and cyclic fatigue, and aging resistance is assessed. We also report on the load to failure of femoral heads prototypes during compression tests. This over-

all reliability assessment ensures a potential future development for these kinds of new

^{*} Corresponding address: Laboratorio di Tecnologia Medica, Istituto Ortopedico Rizzoli, Via di Barbiano, 1/10, 40136 Bologna, Italy. Tel.: +39 051 6366864; fax: +39 051 6366863.

E-mail address: affatato@tecno.ior.it (S. Affatato).

^{1751-6161/\$ -} see front matter C 2010 Elsevier Ltd. All rights reserved. doi:10.1016/j.jmbbm.2010.10.010